

Portable Medium Exchange System



Features

- Automates the process of medium exchange at a specific time interval.
- Battery-operated. Usable in an incubator.
- A standard 6-well plate can be used with this system.
- Wells are observable on a microscope.
- Can be remodeled into a perfusion culture system.

Specifications

Model Number	CEME-0102
Volume of Medium Exchange	Approx. 2.7 ml /cycle
Power Supply	AA (LR6) battery x 2 (for 3 days ~ 7 days)
Wetted Materials	PP, Silicone, Stainless Steel 304
Dimensions	194 x 228 x 140mm

Note: Details including specifications may change without notification.

TAKASAGO ELECTRIC, INC.

66 KAKITSUBATA, NARUMI-CHO, MIDORI-KU, NAGOYA, 458-8522 JAPAN

Tel +81-52-891-2301 Fax+81-52-891-7386

E-mail: info@takasago-elec.co.jp URL: <http://www.takasago-fluidics.com/>